


## Assisted Hearing with the Listening Loop

- The Loop system works *only* with Hearing Aids equipped with a Telecoil.
- If your Hearing aid is controlled by an app be sure to check that the Telecoil function is enabled. This may also be called the T-Coil or Telephone mode.
- Some hearing aids do not have a Telecoil option – therefore the Hearing loop *will not* work with them.
- The Theatre has Headsets that work with the Loop system available from Audience Services in the lobby.



COWICHAN  
**PERFORMING**  
**ARTSCENTRE**

[www.cowichanpac.ca](http://www.cowichanpac.ca)

250-746-2722

