

Connect to the Advanced Listening Hearing Loop

For the induction loop system to function properly, patrons must have a hearing aid or cochlear implant that has a T-Coil (or Telecoil) and the T-Coil is active.



The hearing device must be SWITCHED into Telephone (or “T”) mode. Please contact your audiologist or Hearing Aid practitioner if you are unsure if your device has a T-Coil or instructions on how to change your settings on your personal device.

Once T-mode is enabled, the sounds should feed through the induction system.

Alternative hearing equipment is available should there be any issues connecting to the induction loop. Please ask a staff member if you require assistance.